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Before operating the unit, please read this manual thoroughly and retain it for future reference.

For the customers in the U.S.A.

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

You are cautioned that any changes or modifications not expressly approved in this manual could void your authority to operate this equipment.

All interface cables used to connect peripherals must be shielded in order to comply with the limits for a digital device pursuant to Subpart B of part 15 of FCC Rules.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

For the customers in Canada

CAN ICES-3 (A)/NMB-3(A)

Important safeguards/notices for use in the medical environments

1. All the equipments connected to this unit shall be certified according to Standard IEC60601-1, IEC60950-1, IEC60065 or other IEC/ISO Standards applicable to the equipments.
2. Furthermore all configurations shall comply with the system standard IEC60601-1-1. Everybody who connects additional equipment to the signal input part or signal output part configures a medical system, and is therefore, responsible that the system complies with the requirements of the system standard IEC60601-1-1. If in doubt, consult the qualified service personnel.

3. The leakage current could increase when connected to other equipment.

4. For this particular equipment, all accessory equipment connected as noted above, must be connected to mains via an additional isolation transformer conforming with the construction requirements of IEC60601-1 and providing at least Basic Insulation.

5. This equipment generates, uses, and can radiate radio frequency energy. If it is not installed and used in accordance with the instruction manual, it may cause interference to other equipment. If this unit causes interference (which can be determined by unplugging the power cord from the unit), try these measures: Relocate the unit with respect to the susceptible equipment. Plug this unit and the susceptible equipment into different branch circuit.

Consult your dealer. (According to standard EN60601-1-2 and CISPR11, Class B, Group 1)

This symbol indicates the manufacturer, and appears next to the manufacturer's name and address.

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For the customers in the U.S.A.
SONY LIMITED WARRANTY -

Please visit

<http://www.sony.com/psa/warranty> for important information and complete terms and conditions of Sony's limited warranty applicable to this product.

For the customers in Canada

SONY LIMITED WARRANTY -

Please visit <http://www.sonybiz.ca/pro/lang/en/ca/article/resources-warranty-product-registration> for important information and complete terms and conditions of Sony's limited warranty applicable to this product.

For the customers in Europe

Sony Professional Solutions Europe - Standard Warranty and Exceptions on Standard Warranty. Please visit <http://www.pro.sony.eu/warranty> for important information and complete terms and conditions.

For the customers in Korea

SONY LIMITED WARRANTY -

Please visit

<http://bpeng.sony.co.kr/handler/BPAS-Start> for important information and complete terms and conditions of Sony's limited warranty applicable to this product.

Precaution

On simultaneous use with a radio knife, etc.

If this unit is used together with a radio knife, etc., the picture may be disturbed, warped or otherwise abnormal as a result of strong radio waves or voltages from the device. This is not a malfunction.

When you use this unit simultaneously with a device from which strong radio waves or voltages are emitted, confirm the effect of this before using such devices, and install this unit in a way that minimizes the effect of radio wave interference.

On moisture condensation

If the unit is suddenly taken from a cold to a warm location, or if ambient temperature suddenly rises, moisture may form on the outer surface of the unit and/or inside of the unit. This is known as condensation. If condensation occurs, turn off the unit and wait until the condensation clears before operating the unit. Operating the unit while condensation is present may damage the unit.

On parts life expectancy

The life expectancy of the electrolytic capacitor is about 5 years under normal operating temperatures and normal usage (8 hours per day; 25 days per month). If usage exceeds the above normal usage frequency, the life expectancy may be reduced correspondingly.

GB

Overview

The BKM-250TGM 3G/HD/SD-SDI Input Adaptor is a video signal input adaptor for Sony monitors.

When installed in an input option slot on the main unit, it provides video input and output connectors for the main unit and a decoder for 3G/HD/SD serial digital component signals.

Notes

- There are some monitors that are not compatible with this unit.
To check if the monitor is compatible with this unit, refer to the instructions for use for the monitor.
- The model name displayed on the monitor's STATUS menu is "BKM-250TG."

Function

Decoder for serial digital component signals

The BKM-250TGM supports 3G-SDI that transmits 1080/50P or 1080/60P signals using a single BNC cable, and is equipped with a built-in decoder for 3G/HD/SD serial digital component signal.

Serial digital input and output signal connector


The BKM-250TGM is equipped with two input connectors and monitor output connectors that correspond to each input signal for serial digital signals.

Using the Input and Output Connectors

For information about installing the BKM-250TGM in an input option slot on the main unit, refer to the instructions for use for the monitor.

Configuration of input/output connectors and signals that may be input

The configuration of the input and output connectors, and the signals that may be input are shown below.

Input of serial digital component signals

You can input serial digital signals to connectors ① and output signals displayed on the screen of the main unit or the monitor from connector ②.

You need not attach the 75Ω termination to connector ②.

Notes

- The OUT connector (monitor out) outputs signal only when the power of the main unit is ON. The monitor out signal is not output when the monitor is in standby mode.
- The signal from the OUT connector does not comply with ON-LINE output signal specifications.

Specifications

General

Voltage	+3.3 V, +5 V (supplied from the main unit)
Power consumption	Approx. 4 W
Operating conditions	
Temperature	0 °C to 35 °C (32 °F to 95 °F)
Optimum temperature	20 °C to 30 °C (68 °F to 86 °F)
Humidity	0% to 90% (no condensation allowed)
Pressure	700 hPa to 1060 hPa
Storage and transport temperature	-20 °C to +60 °C (-4 °F to +140 °F)
Storage and transport humidity	0% to 90% (no condensation allowed)
Storage and transport pressure	700 hPa to 1060 hPa
Maximum external dimensions (w/h/d)	100 × 20 × 162 mm (4 × ¹³ / ₁₆ × 6 ¹ / ₂ inches)
Mass	Approx. 270 g (9.5 oz)

Input/Output connectors

Digital input	BNC (2), with monitor output connector
---------------	--

Available signal format

Refer to the instructions for use of the monitor.

Signal characteristics

Digital component signals	
Sampling frequency	
3G-SDI: Y/Cb/Cr (4:2:2)	148.5 MHz/74.25 MHz/ 74.25 MHz
Y/Cb/Cr (4:4:4)	148.5 MHz/148.5 MHz/ 148.5 MHz
G/B/R	148.5 MHz/ 148.5 MHz/148.5 MHz
HD-SDI: Y/Cb/Cr	74.25 MHz/37.125 MHz/ 37.125 MHz
SD-SDI: Y/Cb/Cr	

13.5 MHz/6.75 MHz/
6.75 MHz

Quantization	3G-SDI: 10 bits/sample, 12 bits/ sample
	HD-SDI: 10 bits/sample
	SD-SDI: 10 bits/sample
OUT connector (monitor out)	
Output signal amplitude:	800 mVp-p ±10%
Output impedance:	75Ω unbalanced
Transmission distance	
3G-SDI:	70 m (approx. 230 ft) max. (When using 5C-FB coaxial cables (Fujikura. Inc.) or equivalent.)
HD-SDI:	100 m (approx. 328 ft) max. (When using 5C-FB coaxial cables (Fujikura. Inc.) or equivalent.)
SD-SDI:	200 m (approx. 656 ft) max. (When using 5C-2V coaxial cables (Fujikura. Inc.) or equivalent.)

Note

Not available for analog signal.

Supplied accessories

Instructions for Use (1)
European Representative (1)

Design and specifications are subject to change without notice.

Notes

- Always verify that the unit is operating properly before use. SONY WILL NOT BE LIABLE FOR DAMAGES OF ANY KIND INCLUDING, BUT NOT LIMITED TO, COMPENSATION OR REIMBURSEMENT ON ACCOUNT OF THE LOSS OF PRESENT OR PROSPECTIVE PROFITS DUE TO FAILURE OF THIS UNIT, EITHER DURING THE WARRANTY PERIOD OR AFTER EXPIRATION OF THE WARRANTY, OR FOR ANY OTHER REASON WHATSOEVER.
- SONY WILL NOT BE LIABLE FOR CLAIMS OF ANY KIND MADE BY USERS OF THIS UNIT OR MADE BY THIRD PARTIES.
- SONY WILL NOT BE LIABLE FOR THE TERMINATION OR DISCONTINUATION OF ANY SERVICES RELATED TO THIS UNIT THAT MAY RESULT DUE TO CIRCUMSTANCES OF ANY KIND.

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© 2014 Sony Corporation
(Fujikura. Inc.) of equivalent.)
coaxkabels (Fujikura. Inc.) of